



NATIONAL JUNIOR COLLEGE
SH 2 Year – End Practical Examination
Higher 2

CANDIDATE
NAME

SUBJECT
CLASS

REGISTRATION
NUMBER

CHEMISTRY

9729/04

Paper 4 Practical
Candidate answer on the Question paper.


Tuesday 20 August 2019
2 hours 30 minutes

READ THE INSTRUCTIONS FIRST Write your identification number and name. Circle the practical shift and laboratory where appropriate, in the boxes provided. Write in dark blue or black pen. You may use an HB pencil for any diagrams or graphs. Do not use staples, paper clips, highlighters, glue or correction fluid. Answer all questions in the spaces provided on the Question Paper. The use of an approved scientific calculator is expected, where appropriate. You may lose marks if you do not show your working or if you do not use appropriate units. Qualitative Analysis Notes are printed on pages 19 and 20. At the end of the examination, submit the question paper. The number of marks is given in brackets [] at the end of each question or part question.	Shift			
	1	2	3	
	Laboratory			
	CM41	CM42	CM43	CM44
	PH31	PH32	PH33	
	BI23	BI24		
	For Examiner's Use			
	1	/ 12		
	2	/ 21		
	3	/ 14		
4	/ 8			
Total	/ 55			

This document consists of **19** printed pages and **1** blank page.

1 Determination of the solubility product of calcium iodate(V)

When solid calcium iodate(V) is added to water to form a saturated solution, an equilibrium between the undissolved salt and its aqueous ions is established.




The solubility product, K_{sp} , of calcium iodate(V) is expressed as follows.

$$K_{\text{sp}} = [\text{Ca}^{2+}(\text{aq})] [\text{IO}_3^-(\text{aq})]^2$$

The solubility product can be found by determining the equilibrium concentration of IO_3^- ions in a saturated solution of calcium iodate(V) through a redox titration.

Iodate(V) ions, IO_3^- , react with iodide ions according to the following equation.



The iodine produced in this reaction may be titrated against sodium thiosulfate using starch solution as indicator. Thiosulfate ions react with iodine according to the following equation:



In this question, you will prepare a saturated solution of calcium iodate(V) to carry out the redox titration. Using the experimental results, you will determine the solubility product of calcium iodate(V).

You are provided with:

FA 1 is solid calcium iodate(V), $\text{Ca}(\text{IO}_3)_2$

FA 2 is 0.05 mol dm^{-3} sodium thiosulfate, $\text{Na}_2\text{S}_2\text{O}_3$

FA 3 is 0.50 mol dm^{-3} potassium iodide, KI

FA 4 is 0.50 mol dm^{-3} sulfuric acid, H_2SO_4

FA 5 is starch indicator

Preparation of a saturated solution of calcium iodate (V), FA 6

(a) Procedure

1. Use a 100 cm^3 measuring cylinder to transfer 80.0 cm^3 of deionised water into a 250 cm^3 beaker.
2. Transfer all **FA 1** from the weighing bottle into the beaker. Stir for 1 to 2 minutes and leave to stand for **5 minutes**. There will be some $\text{Ca}(\text{IO}_3)_2$ solids left undissolved.
3. Use a 0.2°C division thermometer, record the temperature of the solution in the beaker.
4. To remove the undissolved solids, filter the saturated solution into another clean and **dry** 250 cm^3 beaker using a **dry** filter paper and **dry** funnel. Label this solution as **FA 6**.

Titration

(b) Procedure

- (i)
1. Fill a burette with **FA 2**.
 2. Pipette 10.0 cm³ of **FA 6** into a conical flask.
 3. Use a 10 cm³ measuring cylinder to add 5.0 cm³ of **FA 3**, followed by 5.0 cm³ of **FA 4** to the same conical flask and swirl quickly.
 4. Titrate the iodine in the conical flask with **FA 2** until the brown colour of the iodine becomes pale yellow.
 5. Add about 10 drops of starch solution to the flask and continue adding **FA 2** until the blue-black colour just disappears.
 6. Repeat the titration as many times as you think necessary to obtain consistent results. Record your titration results in the space below.

Results

Temperature of the solution in the beaker = °C

	1	2	
Final burette reading / cm ³			
Initial burette reading / cm ³			
Volume of FA 2 used / cm ³			

- [1] Record the temperature of the solution in the beaker
 [1] Records all burette readings to 2 d.p including a(ii). (MMO)
 [1] Has at least 2 titre values ± 0.10 cm³ (MMO)
 [1] Table has appropriate headings and units (PDO)

[-1] Incorrect calculation for volume of FA2

[4]

- (ii) From your titres, obtain a suitable volume of **FA 2** to be used in your calculations. Show clearly how you obtained this volume.

Average vol of **FA 2** used = (x.00 + .00)/2 = x.00 cm³ (round to 2 d.p.) [1]

- * Do not award marks if inconsistent values used
 * Do penalise if no working

10.0 cm³ of **FA 6** produced an amount of iodine which required cm³ of **FA 2**.
 [1]

(c) Calculations

Show your working and appropriate significant figures in the final answer to **each** part.

- (i) Calculate the number of moles of $\text{S}_2\text{O}_3^{2-}$ in the volume of **FA 2** recorded in **b(ii)**, and hence determine the molar concentration of IO_3^- ions in **FA 6**.

$$\text{Amt of } \text{S}_2\text{O}_3^{2-} = (\text{ans from b(ii)}/1000) \times 0.05 = A \text{ mol [1]}$$

$$\frac{\eta_{\text{IO}_3^-}}{\eta_{\text{S}_2\text{O}_3^{2-}}} = \frac{1}{6}$$

$$\text{Amt of } \text{IO}_3^- = \frac{1}{6} \times A = B \text{ mol}$$

$$[\text{IO}_3^-] = B \times \frac{1000}{10.0} = c \text{ mol dm}^{-3} \text{ [1] (Allow ecf)}$$

$$\text{number of moles of } \text{S}_2\text{O}_3^{2-} = \dots\dots\dots$$

$$\text{concentration of } \text{IO}_3^- \text{ ions in FA 6} = \dots\dots\dots [2]$$

- (ii) Determine the molar concentration of Ca^{2+} in **FA 6** and hence, calculate the solubility product, K_{sp} , of calcium iodate(V), stating the units.

$$\frac{\eta_{\text{Ca}^{2+}}}{\eta_{\text{IO}_3^-}} = \frac{1}{2}$$

$$[\text{Ca}^{2+}] = \frac{1}{2} [\text{IO}_3^-] = \frac{1}{2} c \text{ mol dm}^{-3} \text{ [1] (Allow ecf)}$$

$$\begin{aligned} K_{\text{sp}} &= [\text{Ca}^{2+}(\text{aq})] [\text{IO}_3^-(\text{aq})]^2 \\ &= [\frac{1}{2} C] [C]^2 \\ &= \frac{1}{2} C^3 \text{ mol}^3 \text{ dm}^{-9} \text{ [1] (Allow ecf)} \\ &\quad \text{(Do not award the mark if the units is wrong)} \end{aligned}$$

$$\text{concentration of } \text{Ca}^{2+} \text{ ions in FA 6} = \dots\dots\dots$$

$$\text{solubility product of } \text{Ca}(\text{IO}_3)_2 = \dots\dots\dots [2]$$

- (d) Explain why it was necessary to use the **dry** apparatus (e.g. beaker, filter paper and funnel) in the preparation of **FA 6**, and what effect failing to do it would have on the titre values?

To prevent the dilution of the saturated **FA 6** solution which results in lesser amount of

IO_3^- reacted during titration and hence smaller titre value [1]

(Do not award mark if there is no mention to the change in the titre value or change in conc of **FA 6** or IO_3^- .)


[1]

- (e) A literature value for the solubility product of calcium iodate(V) is found to be 6.15×10^{-6} at 20 °C.

State a possible reason for the difference in the K_{sp} value that you have calculated and suggest an improvement that might allow a value closer to the literature value to be obtained.

If the calculated K_{sp} > literature K_{sp} value, the explanation is as follow:

The experiment is not carried out at 20 °C [1].



Hence, the equilibrium position shifted to the right.

Use a thermostatically controlled water bath maintained at 20 °C during the preparation of FA 6 before carrying out the filtration [1].

or

If the calculated K_{sp} < literature K_{sp} value, the explanation is as follow:

The salt has not dissolve to form a saturated solution or reach equilibrium with its dissolved ions. [1]

Leave the mixture to stand for longer period and use a thermostatically controlled water bath to maintain at 20 °C for the salt solution to reach equilibrium / saturation. [1]

[Total: 12]

2 Determination of the enthalpy change of neutralisation, ΔH_n

FA 7 is a solution of HCl of unknown concentration

FA 8 is a solution of 1.0 mol dm^{-3} of NaOH

In this question, you will carry out a series of experiments where different volumes of **FA 7** and **FA 8** are mixed together.

You will determine the temperature change of the reaction mixture, ΔT , of each experiment and then analyse your results graphically in order to determine the

- concentration of HCl in **FA 7**
- maximum temperature change, ΔT_{max}
- value for the enthalpy change of neutralisation, ΔH_n

(a) Determining the change in temperature for a series of reactions between **FA 7** and **FA 8**

(i) Experiment 1

1. Place a Styrofoam cup inside a second Styrofoam cup held in a 250 cm^3 beaker to prevent it from tipping over.
2. Using a measuring cylinder, transfer 20.0 cm^3 of **FA 7** into the Styrofoam cup.
3. Measure the temperature of **FA 7** in the Styrofoam cup. Record the initial temperature of the solution of **FA 7** as $T_{\text{FA 7}}$.
4. Rinse and dry the thermometer.
5. Measure 50.0 cm^3 of **FA 8** using another measuring cylinder. Record the initial temperature of the solution of **FA 8** as $T_{\text{FA 8}}$.
6. Transfer the **FA 8** into the cup containing **FA 7**. Stir and record the maximum temperature of the reaction as T_{max} .
7. Rinse and dry the Styrofoam cup and thermometer.

Experiment 2

Repeat steps 2 to 7 using 50.0 cm^3 of **FA 7** and 20.0 cm^3 **FA 8**.

Experiments 3 to 6

Carry out **four** further experiments to investigate how the T_{max} changes with different volumes of **FA 7** and **FA 8**. In each case, the total volume of the reaction mixture must be 70.0 cm^3 and the volume of each reagent should not be less than 20.0 cm^3 .

In an appropriate format in the space provided on next page, record

- all measurements of volumes used,
- all temperatures measured, $T_{\text{weighted initial}}$ and ΔT

$$\Delta T = T_{\text{max}} - T_{\text{weighted initial}}$$

$$T_{\text{weighted initial}} = \frac{(\text{Volume of FA 7} \times T_{\text{FA 7}}) + (\text{Volume of FA 8} \times T_{\text{FA 8}})}{\text{Volume of FA 7} + \text{Volume of FA 8}}$$

Result

Expt	Vol of FA 7 /cm ³	Vol of FA 8 /cm ³	T _{FA 7} / °C	T _{FA 8} / °C	T _{weighed initial} / °C	T _{may} / °C	ΔT / °C
1							
2							
3							
4							
5							
6							

Results							
Expt	Volume of FA7/cm ³	T _{FA7} / °C	Volume of FA8/cm ³	T _{FA8} / °C	T _{weighed initial} / °C	T _{max} / °C	ΔT / °C
1	20.0	29.4	50.0	30.0	29.8	32.8	3.0
2	50.0	29.8	20.0	30.0	29.9	32.6	2.7
3	30.0	29.8	40.0	30.0	29.9	34.4	4.5
4	40.0	29.8	30.0	29.9	29.8	35.3	5.5
5	35.0	29.8	35.0	30.0	29.9	35.1	5.2
6	45.0	29.8	25.0	28.9	29.8	34.4	4.6

[1] Appropriate header with units

[1] Recording to correct precision

[1] Correct calculation of T_{weighed initial} and ΔT_{max} in 1 d.p

[1] Proper spread of the choice of the volume used for the other 4 expts

[-1] Missing any of the experiment results

[-1] Missing the recording of T_{weighed initial}[-1] Vol of either FA7 or FA8 used is less than 20.0 cm³[-1] Recording of data is NOT IN ONE table.[-1] Missing T_{FA7} and/or T_{FA8}

[4]

- (ii) Plot a graph of ΔT (y-axis) against volume of **FA 7** used (x-axis) on the grid in Fig. 2.1 using the data you have obtained in **2a(i)**.




Fig. 2.1

[3]

[1] - Label of the axis with units and appropriate scale

[1] - Points plotted within half a square of 'x'

[1] - The line of best fit must include two points at $x = 0$ and $x = 70$

Note: $\Delta T_{\max}$ should be higher than any of the plotted points.

[-1] if $\Delta T_{\max}$ is lower than the plotted points.

Plotted points covers more than 50% of grids in both x & y axis

- (iii) Draw two best-fit straight lines, the first should be drawn using the plotted points where ΔT is increasing and the second should be drawn using the plotted points where ΔT is decreasing. Extrapolate (extend) both lines until they intercept.

[1]

[1] - Line of best fit to pass through the plotted points

- (iv) From your graph in Fig. 2.1, determine the maximum temperature change, $\Delta T_{\max}$, and the volume, $V_{\max}$, of **FA 7** required to obtain this.

Correct reading off the $\Delta T_{\max}$ value from the graph [1] (ignore precision)

Correct reading off the $V_{\max}$ value from the graph [1] (ignore precision)

$$\Delta T_{\max} = \dots\dots\dots^{\circ}\text{C} \quad V_{\max} = \dots\dots\dots\text{cm}^3 \quad [2]$$

(b) Calculation

Show your working and appropriate significant figures in the final answer to **each** part.

Using your answers in **2a(iv)**, calculate

(i) the concentration, in mol dm^{-3} , of HCl in **FA 7**.

$$V_{\text{NaOH}} = 70 - V_{\max} = A \text{ cm}^3$$

$$\text{Amt of NaOH} = A/1000 \times 1.0 = B \text{ mol [1]}$$

$$\text{Amount of HCl} = B \text{ mol}$$

$$\text{Conc of HCl} = B / (V_{\max}/1000) = C \text{ mol dm}^{-3} \text{ [1] (Allow ecf)}$$

$$\text{concentration of HCl in FA 7} = \dots\dots\dots [2]$$

(ii) the heat change for the neutralisation at $\Delta T_{\max}$.

You may assume that the specific heat capacity of the reaction mixture is $4.18 \text{ J g}^{-1} \text{ K}^{-1}$ and that the density is 1.00 g cm^{-3} .

$$q = mc\Delta T_{\max} = 70 \times 4.18 \times \Delta T_{\max} = D \text{ J}$$

$$\text{heat change} = -D \text{ J [1] (Accept if heat change value is positive)}$$

$$\text{heat change} = \dots\dots\dots [1]$$

(iii) Hence, calculate the enthalpy change of neutralisation, ΔH_n .

$$\text{Amount of H}_2\text{O} = \text{Amount of NaOH} = \text{Amount of HCl} = B \text{ mol}$$

$$\Delta H_{\text{neut}} = -D / B = -E \text{ kJ mol}^{-1} \text{ [1] (Must have the negative sign to award the mark) (Allow ecf)}$$

$$\Delta H_n = \dots\dots\dots [1]$$

(c) In another similar experiment, 1.0 mol dm^{-3} aqueous ammonia was used instead of 1.0 mol dm^{-3} NaOH. A graph of maximum change in temperature, ΔT , against volume of **FA 7** was plotted.

Suggest how would the $\Delta T_{\max}$ and volume, $V_{\max}$, of **FA 7** from this experiment be compared to your result obtained in **2a(iv)**.

Lower temperature rise, $\Delta T_{\max}$, as some heat released is used to dissociate the weak
base completely. [1]


(Do not accept rate of reaction is less vigorous and heat lost or less water formed or the reaction is less exothermic)

.....
 $V_{\max}$ is the same as same number of moles of HCl is needed to react with the aqueous


 NH_3 since the reacting mole ratio remains the same. [1]

.....[2]

- (d) A student employed another method to determine the concentration of the HCl in **FA 7** through the use of conductometry. Conductometry is the measurement of the electrolytic conductivity to monitor the progress of a chemical reaction due to ions being produced or reacted away. The equivalence point can then be determined for the acid – base reaction by plotting the conductivity against the volume of the HCl added.



The following apparatus were set up.



The student measured the conductivity of the solution in the beaker with each addition of a small fixed volume of HCl and obtained the following graph.





Fig. 2.2

- (i) Indicate on the graph (Fig. 2.2), the point which marks the volume of HCl needed for complete reaction. [1]

$V_{\max}$ is the volume of HCl that corresponds to the minimum conductivity [1]



- (ii) Explain the shape of the graph in Fig. 2.2.

The conductivity decreases as more HCl is added. This is due to the **formation of water** from H^+ and OH^- , which **decreases the total amount of ions** in water and hence conductivity. [1]

The conductivity decreases to the lowest point where the equivalence volume of HCl is reached, as **all the NaOH has reacted** with the HCl. Thus, that is the **minimum amount of ions** that can be present and gives the lowest conductivity. [1]

After the equivalence point, **no more neutralisation** occurs. Any HCl that is added **dissociates and contributes to the total amount of ions** in water, thus increasing the conductivity. [1]

[3]

- (iii) With reference to your experiment in 2(a), suggest which method (calorimetry or conductometry) is more accurate for the determination of the concentration of HCl in FA 7.

Conductometry is more accurate. For conductometry, there is no heat loss to the surroundings.

Explanation 1: For calorimetry, the heat loss to the surroundings will lower the $\Delta T_{\max}$ and hence render the volume of HCl needed for complete neutralisation inaccurate

Explanation 2: For calorimetry, the heat loss to the surroundings will affect each ΔT plotted point, resulting in inaccurate two best-fit straight lines drawn to determine $\Delta T_{\max}$ and hence volume of HCl needed for complete neutralisation will be inaccurate.

[1]

[1]

[Total: 21]

3 Identification of ions

FA 9 is a solution containing a mixture of two salts.

You are to perform the tests described in Table 3.1 and record your observations in the table.

In all tests, the reagents should be added gradually until no further change is observed unless you are instructed otherwise. You should indicate clearly at which stage in a test a change occurs, recording your observations alongside the relevant tests.

The volumes given below are approximate and should be estimated rather than measured, unless otherwise stated. If there is no observable change, write **no observable change**.

No additional or confirmatory tests for ions present should be attempted.

Table 3.1

	<i>Test</i>	<i>Observation</i>
(a)	Place about 3 cm depth of FA 9 in a test-tube, add equal volume of aqueous sodium hydroxide. Warm the contents in the tube.	White, off-white, buff or light brown ppt. [1] Moist red litmus did not turn blue or no alkaline gas produced. [1] Ppt turned browned or darkened. [1]
(b)	Filter the mixture from (a) and collect the filtrate in a test tube. Label the filtrate as FA 10 and retain it for tests 3(c) and 3(d). Leave the residue in the filter paper and observe it again after 5 minutes.	White, off-white, buff, light brown or brown residue. [awarded once either here or in test (a) but not in both] Colourless filtrate obtained. [1] Residue turned brown or darkened. [awarded once either here or in test (a) but not in both]
(c)	Place 2 cm depth of FA 10 in a test-tube and add nitric acid until no further change is seen. Label the solution as FA 11 and retain it for test 3(e).	White ppt formed [1] Dissolved in excess nitric acid to give a colourless solution. [1]
(d)	Place 1 cm depth of FA 10 in a test-tube. Add half a spatula of zinc powder. Warm the contents in the tube cautiously.	Colourless and pungent gas evolved [1] which turned moist red litmus paper blue. [1]
(e)	Place 2 cm depth of FA 11 in a test-tube, add equal depth of aqueous ammonia. <u>Add</u> excess dilute nitric acid followed by aqueous silver nitrate.	White ppt formed [1] Soluble/dissolved in excess aqueous ammonia to form a colourless solution. [1] On adding HNO ₃ acid, no gas evolved [1] White ppt formed [1] Dissolved in excess nitric acid to give a colourless solution. [1] On adding silver nitrate, no ppt observed. [1]

	Test	Observation
(f)	Place 1 cm depth of FA 9 in a test-tube. Add dilute hydrochloric acid followed by aqueous barium chloride.	On adding acid, no observable change/ no gas liberated/ no ppt formed. [1] On adding barium chloride, white ppt formed. [1]

[6]

[6] – obtain 13 points

[5] – obtain 11 points

[4] – obtain 9 points

[3] – obtain 7 points

[2] – obtain 5 points

[1] – obtain 3 points

- (g) (i) Identify the
- two**
- cations and
- two**
- anions present in
- FA 9**
- .

 Mn^{2+} [1] Zn^{2+} [1]

cation 1:

cation 2:



 SO_4^{2-} [1] NO_3^- [1]

anion 1:

anion 2:

[4]

- (ii) Write equations, with state symbols, to illustrate the observations in 3(c).



.....[2]

- (iii) Predict what will happen when solid ammonium nitrate is added to
- FA 9**
- followed by aqueous ammonia.



The dissociation of ammonium nitrate produces NH_4^+ which suppressed the dissociation of NH_3 (or causes the position of equilibrium of (1) lies to the left). [1]

The concentration of OH^- is too low for the ppt to form or less ppt formed. [1]

.....[2]


[Total : 14]

4 Planning

Determination of the hardness of water in a sample.

Hard water is water that has high mineral content and is formed when water percolates through deposits of limestone, chalk and gypsum, which are largely made of calcium carbonates, bicarbonates and sulfates. The hardness of water is measured in terms of the molar concentration of Ca^{2+} in mol dm^{-3} .

In a particular sample of water, the hardness is caused by the presence of calcium chloride. The amount of Ca^{2+} , and hence the concentration of Ca^{2+} , can be determined by precipitation with an excess of oxalic acid.



You are to design an experiment to determine the hardness of water in a sample.

You are provided with

- sample of water with estimated hardness of 0.02 mol dm^{-3}
 - 0.1 mol dm^{-3} of oxalic acid $\text{H}_2\text{C}_2\text{O}_4$
 - electronic balance
 - common laboratory apparatus
- (a) Given the K_{sp} of calcium oxalate is $2.7 \times 10^{-9} \text{ mol}^2 \text{ dm}^{-6}$, predict what will be observed when 25.0 cm^3 of the water sample is mixed with equal volume of the 0.1 mol dm^{-3} of oxalic acid solution.

$$\begin{aligned} [\text{Ca}^{2+}] &= \frac{1}{2}(0.02) = 0.01 \text{ mol dm}^{-3} \\ [\text{C}_2\text{O}_4^{2-}] &= \frac{1}{2}(0.1) = 0.05 \text{ mol dm}^{-3} \\ \text{ionic product, } IP &= [\text{Ca}^{2+}][\text{C}_2\text{O}_4^{2-}] \\ &= (0.01)(0.05) \\ &= 5.0 \times 10^{-4} \text{ mol}^2 \text{ dm}^{-6} [1] \end{aligned}$$

Since $IP > K_{\text{sp}}$, calcium oxalate will be precipitated. [1] (Can award this mark based on correct deduction of IP against K_{sp} even if IP above is calculated wrongly)

[2]

- (b) Your plan should include a step-by-step description of the method, including the
- sequence of steps to ensure that the maximum amount of precipitate is recovered,
 - the use of appropriate apparatus,
 - the recordings using the letters A, B, C etc to represent the measurement taken and present it with an appropriate table format,
 - the processing of the results to determine the molar concentration of Ca^{2+} in the water sample.

[A_r: C, 12.0; Ca, 40.1; H, 1.0; O, 16.0]

Forming of precipitate

1. Pipette out 25.0 cm³ of the sample water with a 25.0 cm³ pipette into a dry and clean 250 cm³ beaker
2. Using a burette, transfer 25.0 cm³ of 0.1 mol dm⁻³ of oxalic acid into the same beaker (as long as oxalic acid is in excess)
3. Stir the mixture thoroughly with a glass rod and leave it to stand for a few minutes for equilibrium to be reached

Recovering of precipitate

4. Filter the solution into a 250 cm³ conical flask and collect the residue precipitate into a pre-weighed filter paper.
5. Rinse the beaker and glass rod with small volume of deionised water and transfer the washings and continue to filter.
6. Wash the precipitate with limited cold deionised water to get rid of aq ions adsorbed on the surface of precipitate.
7. Drain off excess water by pressing the solid between filter papers.
8. Dry the precipitate in oven at 80 - 100 °C $\pm$ 10%.
9. Cool and weigh the dried precipitate
10. Repeat the heat cool weigh process until mass is constant to $\pm 0.01\text{g}$

Mass of filter paper	/ g	A
Mass of filter paper and residue after first heating	/ g	B
Mass of filter paper and residue after second heating	/ g	C
Mass of filter paper and residue after third heating	/ g	C
Mass of residue	/g	C – A = E

$$\text{Amount of CaC}_2\text{O}_4 = \frac{E}{40.1 + 2(12.0) + 4(16.0)} = F \text{ mol}$$



$$\text{Amt of Ca}^{2+} \text{ to form CaC}_2\text{O}_4 = F \text{ mol}$$

$$[\text{Ca}^{2+}] \text{ in mol dm}^{-3} = \frac{F}{\left(\frac{25.0}{1000}\right)} = G \text{ mol dm}^{-3}$$

Marks	Requirement	Key marking Points
L-[1]	Logical flow of plan	- Pipetting out water sample into beaker followed by adding fixed solution of oxalic acid into beaker - Filtering out precipitate - Drying with filter paper and oven or IR lamp (Buchner funnel accepted) <i>Award mark for including all 3 points.</i>
A- [1]	Appropriate apparatus and their <u>capacities</u>	100 cm³ beaker 250 cm³ beaker 250 cm³ conical flask 25.0 cm³ pipette
R-[1]	Essential <u>experimental details</u> to ensure reliable results	- Stir mixture with glass rod and leave to stand for a few minutes - Repeat heat cool weigh process until mass is constant at $\pm 0.01\text{g}$
P-[1]	Appropriate volumes with correct precision	25.0 cm³ of water sample 25.0 cm³ of oxalic acid
D-[1]	Appropriate table for tabulation of data	- Correct headers and units - Letter
M-[1]	Manipulation of data	- Correct calculation for mass of residue - Correct calculation for the conc of Ca²⁺ in water

[6]

[Total: 8]